

evochrom

Diesel System for EN 12916, IP 391, ASTM D6591

- Complete HPLC Pack system allows plug and go for the busy lab
- Precision engineered HPLC hardware with proven German technology
- Chromatography assured with optimised column technology meeting the latest EN 12916, IP 391 and ASTM D6591 Method requirements
- Quality separations to provide accurate and precise results every time
- Fully compatible to meet and measure to EN 590 Diesel Fuel Specifications for Aromatics

evochrom Diesel System includes:

- Sample Injector
- HPLC Isocratic Pump
- Optional Solvent Reservoir/Degasser Unit
- Column Oven (with Safety Solvent Alarm System)
- Amino Column as specified in Test Method
- Refractive Index Detector
- ChromGate Integration Software Package
- STVM Validation Materials



Diesel Chromatogram
by IP 391, EN 12916 & ASTM D6591

