

evochrom

Jet/Diesel Fuel Static Dissipator Additive (SDA) System for IP PM DP, ASTM Draft

- Complete HPLC Pack system allows plug and go for the busy lab
- Precision engineered HPLC hardware with proven German technology
- Chromatography assured with optimised column technology meeting the latest IP PM DP and ASTM Draft Method requirements
- Quality separations to provide accurate and precise results every time
- Ideal for QC checking in distribution systems

evochrom SDA System includes:

- Sample Injector
- Solid Phase Extraction (SPE) Cartridges
- HPLC Isocratic Pump
- Optional Solvent Reservoir/Degasser Unit
- C6 alkyl boned RP Column as specified in Test Method
- Variable Wavelength UV Detector
- ChromGate Integration Software Package



SDA Chromatogram by IP PM DP & ASTM Draft Method

